



1700 Elm Street Minneapolis, MN 55414

Phone: 612.607.1700 Fax: 612.607.6444

DETERMINATION OF PCDD/PCDF LEVELS

Prepared for: **ESS Laboratory** Attn: Jena Paola 185 Frances Avenue Cranston, RI 02910-2211



This report contains 25 pages.

The results reported herein conform to the most current NELAC standards, where applicable, unless otherwise narrated in the body of the report.


Project: Chemical Analysis

Client Project Number: 0606373

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full, without the written consent of Pace Analytical Services, Inc.





1700 Elm Street Minneapolis, MN 55414

> Phone: 612.607.1700 Fax: 612.607.6444



REPORT OF: CHEMICAL ANALYSES

PROJECT:

PCDD/PCDF ANALYSES

DATE: July 12, 2006

ISSUED TO:

ESS Laboratory

REPORT NO: 06-1034231

Attn: Jena Paola 185 Frances Avenue

Cranston, RI 02910-2211

INTRODUCTION

This report presents the results from the analyses performed on ten samples submitted by a representative of ESS Laboratory. The samples were analyzed for the presence or absence of polychlorodibenzo-p-dioxins (PCDDs) and polychlorodibenzofurans (PCDFs) using a modified version of USEPA Method 8290.

SAMPLE IDENTIFICATION

Client ID	Sample Type	Date Received	PACE ID
0606373-01	Solid	06/23/06	1034231001
0606373-03	Solid	06/23/06	1034231002
0606373-05	Solid	06/23/06	1034231003
0606373-07	Solid	06/23/06	1034231004
0606373-09	Solid	06/23/06	1034231005
0606373-12	Solid	06/23/06	1034231006
0606373-14	Solid	06/23/06	1034231007
0606373-16	Solid	06/23/06	1034231008
0606373-18	Solid	06/23/06	1034231009
0606373-20	Solid	06/23/06	1034231010

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full, without the written consent of Pace Analytical Services, Inc.





Pace Analytical Services, Inc. 1700 Elm Street

Minneapolis, MN 55414

Phone: 612,607,1700 Fax: 612.607.6444

www.pacelabs.com REPORT OF: CHEMICAL ANALYSES

PROJECT:

PCDD/PCDF ANALYSES

DATE: July 12, 2006

PAGE: 2

REPORT NO: 06-1034231

RESULTS

The results are included in the following:

Appendix A – Chain of Custody Documentation

Appendix B - PCDD/PCDF Results

DISCUSSION

The recoveries of the isotopically-labeled PCDD/PCDF internal standards in the sample extracts ranged from 34-142%. With the exception of seven values, which were flagged "P" on the sample results tables; the labeled standard recoveries obtained for the samples were within the 40-135% target range specified in Method 8290. Also, since the quantification of the native 2,3,7,8-substituted congeners was based on isotope dilution, the data were automatically corrected for variation in recovery and accurate values were obtained.

In some cases, interfering substances impacted the determinations of PCDD or PCDF congeners. The affected values were flagged "I" where incorrect isotope ratios were obtained, or "E" where polychlorinated diphenyl ethers were present.

A laboratory method blank was prepared and analyzed with each sample batch as part of our routine quality control procedures. The results show the blanks, with the exception of a trace level of OCDD, to be free of PCDDs and PCDFs at the reporting limits. The OCDD was below the calibration range of the method. One sample contained a similar level to the corresponding blank and was flagged "B" on the results table. This value may be, at least partially, attributed to the background. It should be noted that levels less than ten times the background are not generally considered to be statistically different from the background.

Laboratory and matrix spike samples were also prepared with the samples using clean sand or sample material that had been fortified with native standards. The results show that the spiked native compounds were recovered at 81-108% in the lab spikes and at 57-113% in the background subtracted matrix spikes. Relative percent differences in the matrix spikes ranged from 1.0-11%. These results indicate high degrees of accuracy and precision for these determinations.

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full. without the written consent of Pace Analytical Services, Inc.





Pace Analytical Services, Inc. 1700 Elm Street Minneapolis, MN 55414

> Phone: 612.607.1700 Fax: 612.607.6444

REPORT OF: CHEMICAL ANALYSES

PROJECT:

PCDD/PCDF ANALYSES

DATE: July 12, 2006

PAGE: 3

REPORT NO: 06-1034231

DISCUSSION (Cont.)

The responses obtained for several compounds in ending calibration standards F60629A_16 and F60708A_17 were outside the target range for this method. As specified in the method, the averages of the daily response factors for these compounds were used in the calculations for the samples from these run-shifts. It should be noted that the accuracy of the native congener determinations was not affected by these deviations.

REMARKS

The sample extracts will be retained for a period of 15 days from the date of this report and then discarded unless other arrangements are made. The raw mass spectral data will be archived on magnetic tape for a period of not less than one year. Questions regarding the data contained in this report may be directed to the author at the number provided below.

Pace Analytical Services, Inc.

Scøtt 🗘 Unze

Project Manager, HRMS

(612) 607-6383

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full, without the written consent of Pace Analytical Services, Inc.





Phone: 612.607,1700 Fax: 612.607.6444



TABLE 1. 2,3,7,8-TCDD Equivalency Factors (TEFs) for the Polychlorinated Dibenzo-p-dioxins and Dibenzofurans

Number	Compound(s)	TEF
1	2,3,7,8-TCDD	1.00
2	1,2,3,7,8-PeCDD	0.50
3	1,2,3,6,7,8-HxCDD	0.1
4	1,2,3,7,8,9-HxCDD	0.1
5	1,2,3,4,7,8-HxCDD	0.1
6	1,2,3,4,6,7,8-HpCDD	0.01
7	OCDD	0.001
8	* Total - TCDD	0.0
9	* Total - PeCDD	0.0
10	* Total - HxCDD	0.0
11	* Total - HpCDD	0.0
12	2,3,7,8-TCDF	0.10
13	1,2,3,7,8-PeCDF	0.05
14	2,3,4,7,8-PeCDF	0.5
15	1,2,3,6,7,8-HxCDF	0.1
16	1,2,3,7,8,9-HxCDF	0.1
17	1,2,3,4,7,8-HxCDF	0.1
18	2,3,4,6,7,8-HxCDF	0.1
19	1,2,3,4,6,7,8-HpCDF	0.01
20	1,2,3,4,7,8,9-HpCDF	0.01
21	OCDF	0.001
22	* Total - TCDF	0.0
23	* Total - PeCDF	0.0
24	* Total - HxCDF	0.0
25	* Total - HpCDF	0.0

^{*}Excluding the 2,3,7,8-substituted congeners.

Reference: International Toxic Equivalence

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full, without the written consent of Pace Analytical Services, Inc.





Pace Analytical Services, Inc. 1700 Elm Street Minneapolis, MN 55414

> Phone: 612.607.1700 Fax: 612.607.6444

APPENDIX A

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full, without the written consent of Pace Analytical Services, Inc.



ESS Laboratory Pace

185 Frances Avenue, Cranston, RI 02910-2211 Tel. (401) 461-7181 Fax (401) 461-4486 Division of Thielsch Engineering, In www.esslaboratory.com

CHAIN OF CUSTODY

ESS LAB PROJECT 1D 3606373 103423/ Page_of Other Format: Excel Access PDF Electronic Deliverable Reporting Limits If faster than 5 days, prior approval by laboratory is required #_ Other, Is this project for any of the following: MA-MCP Navy USACE ME State where samples were collected from:

Ϋ́

 \overline{z}

Ξ

IJ

MA RI

Container Type: P-Poly G-Glass S-Sterile V-VOA Matrix: (S-Soil SD-Solid D-Sludge WW-Waste Water GW-Cround Water SW-Surface Water DW-Drinking Water O-Oil W-Wipes F-Filters 103423100 Preservation Code: 1- NP, 2- HCI, 3- H.SO., 4- HNO., 5- NaOH, 6- MeOH, 7- Asorbic Acid, 8- ZnAct, 9-Write Required Analysis Pres Code Project Name (20 Char. or less) Sample Identification (20 Char. or less) 107 Email Address PO# ĺ Zip Address internal Use Only ر ح MATRIX **JMOC** Fax# Sono Packa ソイかい Collection 2160 780 1015 Yes ____ No 1000 201/6/ 70/2-1 Date Contact Person Cooler Present Telephone # ESS LAB Sample# Co. Name City

Sampled by: [] Pickup No NA: Yes Seals Intact

Received by: (Signature) B. Flow *By circling MA-MCP, client acknowledges samples were collected in accordance with MADEP CAM VII A Finquished by: (Signature)

Please fax all changes to Chain of Custody in writing.

(White) Lab Copy 2 (Yellow) Client Receipt

Date/Time

Received by: (Signature)

Date/Time

Received by: (Signature)

Date/Time

Relinquished by: (Signature)

Date/Time 7, 84 Relinquished by: (Signature)

Received by: (Signature)

Date/Time

Relingurshed by: (Signature)

Cooler Temp:

[] Technicians_

6/2/en 0.15

Date/Time



Pace Analytical Services, Inc. 1700 Elm Street Minneapolis, MN 55414

> Phone: 612.607.1700 Fax: 612.607.6444

APPENDIX B

REPORT OF LABORATORY ANALYSIS

This report shall not be reproduced, except in full, without the written consent of Pace Analytical Services, Inc.





Method 8290 Blank Analysis Results

Client - ESS Laboratory

Lab Sample ID Filename

Total Amount Extracted **ICAL Date** CCal Filename(s)

BLANK-10052 U60702A 03 10.2 g 07/01/2006

U60701A_23 & U60702A_14

Matrix Dilution Extracted

Analyzed

Solid NA 06/27/2006

07/02/2006 10:28

Injected	Bv	BAL

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	ND ND	Mar 100 Mar 100 Mar	0.200 0.200	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	68 63 68
2,3,7,8-TCDD Total TCDD	ND ND	100 mp mp (mp	0.200 0.200	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	73 79 84
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	ND ND ND		0.980 0.980 0.980	1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C 1,2,3,4,7,8-HxCDD-13C	2.00 2.00 2.00 2.00 2.00	74 77 77 77 80
1,2,3,7,8-PeCDD Total PeCDD	ND ND		0.980 0.980	1,2,3,4,7,8-HXCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	63 68 61
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF	ND ND ND ND		0.980 0.980 0.980 0.980	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C 1,2,3,4-TCDD-13C	2.00 4.00 2.00	73 68 NA
Total HxCDF 1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	ND ND ND ND		0.980 0.980 0.980 0.980 0.980	1,2,3,7,8,9-HxCDD-13C 2,3,7,8-TCDD-37Cl4	2.00 0.20	NA 60
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	ND ND ND	100 also 100 also 100 100 also 100 also 100 100 also 100 also 100	0.980 0.980 0.980	Total 2,3,7,8-TCDD Equivalence: 0.00 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	ND ND		0.980 0.980			
OCDF OCDD	ND ND	Mar also had and made	2.000 2.000			

Conc = Concentration (Totals include 2,3,7,8-substituted isomers).

EMPC = Estimated Maximum Possible Concentration

LRL = Lower Reporting Limit

J = Concentration detected is below the calibration range

P = Recovery outside of target range

A = Detection Limit based on signal-to-noise measurement

I = Interference

E = PCDE Interference

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No.....1034231



Method 8290 Blank Analysis Results

Client - ESS Laboratory

Lab Sample ID Filename

Total Amount Extracted ICAL Date

CCal Filename(s)

BLANK-10057 F60630A 04 10.2 g

05/31/2006 F60630A_01 & F60630A_13 Matrix Dilution Extracted

Solid NA 06/28/2006

Analyzed 06/30/2006 11:34 Injected By

CVS

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	ND ND	20 22 42 42 42 42 44 44 44 44 44 44 44 44	0.200 0.200	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	76 76 73
2,3,7,8-TCDD Total TCDD	ND ND		0.200 0.200	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	79 103 75
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	ND ND ND		0.980 0.980 0.980	1,2,3,4,7,6-1;XCDF-13C 1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00 2.00	73 70 79 75
1,2,3,7,8-PeCDD Total PeCDD	ND ND		0.980 0.980	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C	2.00 2.00 2.00	78 76 70
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF	ND ND		0.980 0.980	1,2,3,4,7,8,9-HpCDF-13C 1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 2.00 4.00	59 79 61
2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF Total HxCDF	ND ND ND		0.980 0.980 0.980	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	ND ND ND ND		0.980 0.980 0.980 0.980	2,3,7,8-TCDD-37Cl4	0.20	78
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	ND ND ND		0.980 0.980 0.980	Total 2,3,7,8-TCDD Equivalence: 0.0056 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	ND ND	97 37 97 48 38 97 38 70 48 38	0.980 0.980			
OCDF OCDD	ND 5.6	And 100 440 140 An	2.000 2.000 J			

Conc = Concentration (Totals include 2,3,7,8-substituted isomers).

EMPC = Estimated Maximum Possible Concentration

LRL = Lower Reporting Limit

J = Concentration detected is below the calibration range

P = Recovery outside of target range

A = Detection Limit based on signal-to-noise measurement

I = Interference

E = PCDE Interference

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID Lab Sample ID Filename Injected By **Total Amount Extracted** 0606373-01 1034231001 F60629A_13 BAL

% Moisture

20.9 g 71.2 6.02 g Matrix Dilution

Solid 10

Dry Weight Extracted ICAL Date CCal Filename(s) Method Blank ID

05/31/2006 F60629A 01 & F60629A 16 BLANK-10052

Collected Received Extracted

06/21/2006 06/23/2006 06/27/2006

Analyzed 06/30/2006 00:11

				•		
Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	76 15000		27.0 A 3.3	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C	2.00 2.00	66 70
2,3,7,8-TCDD Total TCDD	33 640		9.7 A 3.3	1,2,3,7,8-PeCDF-13C 2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C	2.00 2.00 2.00	60 68 84
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	230 1600 40000		20.0 AN2 32.0 A 17.0	1,2,3,4,7,8-HxCDF-13C 1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00 2.00	89 45 52 53
1,2,3,7,8-PeCDD Total PeCDD	120 1600		17.0 17.0	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C	2.00 2.00 2.00	63 49 50
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	320 750 640	*****	17.0 17.0 17.0	1,2,3,4,7,8,9-HpCDF-13C 1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 2.00 4.00	39 P 53 38 P
1,2,3,7,8,9-HxCDF Total HxCDF	420 25000	*** ** ** **	17.0 17.0	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	55 150 78 2000		17.0 17.0 17.0 17.0	2,3,7,8-TCDD-37Cl4	0.20	75
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	710 170 2000	~~~~	17.0 17.0 17.0	Total 2,3,7,8-TCDD Equivalence: 1100 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	430 940		17.0 17.0			
OCDF OCDD	190 1600		33.0 33.0			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated * = See Discussion

Report No....1034231



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-03 Lab Sample ID 1034231002 Filename F60629A 09 Injected By BAL **Total Amount Extracted** 20.0 g

Matrix Solid % Moisture 28.7 Dilution NA Dry Weight Extracted 14.3 g Collected 06/21/2006

ICAL Date 05/31/2006 Received 06/23/2006 CCal Filename(s) F60629A 01 & F60629A 16 Extracted 06/27/2006 Method Blank ID BLANK-10052 Analyzed 06/29/2006 20:52

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	4.70 170.00	** *** *** ***	0.54 A 0.14	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	80 87 67
2,3,7,8-TCDD Total TCDD	0.62 21.00		0.53 JA 0.14	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00	78 103
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	76.00 590.00	21 	1.40 EA 0.70 0.70	1,2,3,4,7,6-HXCDF-13C 1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C 1,2,3,4,7,8-HxCDD-13C	2.00 2.00 2.00 2.00 2.00	84 56 62 62 72
1,2,3,7,8-PeCDD Total PeCDD	5.20 55.00		0.70 0.70	1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	58 54 48
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF	18.00 14.00 13.00 9.60		0.70 0.70 0.70 0.70	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C 1,2,3,4-TCDD-13C	2.00 2.00 4.00	70 35 P NA
Total HxCDF	380.00		0.70	1,2,3,7,8,9-HxCDD-13C	2.00	NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	3.40 12.00 7.20 130.00		0.70 J 0.70 0.70 0.70	2,3,7,8-TCDD-37Cl4	0.20	91
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	36.00 6.10 94.00	an and the state of	0.70 0.70 0.70	Total 2,3,7,8-TCDD Equivalence: 51 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	66.00 130.00		0.70 0.70			
OCDF OCDD	37.00 540.00		1.40 1.40			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-05
Lab Sample ID 1034231003
Filename F60629A_14
Injected By BAL
Total Amount Extracted 21.0 g

Total Amount Extracted 21.0 g Matrix Solid % Moisture 73.5 Dilution 10

 Dry Weight Extracted
 5.58 g
 Collected
 06/21/2006

 ICAL Date
 05/31/2006
 Received
 06/23/2006

 CCal Filename(s)
 F60629A_01 & F60629A_16
 Extracted
 06/27/2006

Method Blank ID BLANK-10052 Analyzed 06/30/2006 01:01

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	84 8800	TO 60 TO 10 10	16.0 A 3.6	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00	82 87
2,3,7,8-TCDD Total TCDD	22 420		5.3 A 3.6	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C	2.00 2.00 2.00	69 84 104
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	3100 24000	740	18.0 E 18.0 18.0	1,2,3,4,7,8-HxCDF-13C 1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C 1,2,3,4,7,8-HxCDD-13C	2.00 2.00 2.00 2.00	108 I 60 70 69
1,2,3,7,8-PeCDD Total PeCDD	76 950	160 160 160 160 160 160 160 160	18.0 18.0	1,2,3,4,7,8-HXCDD-13C 1,2,3,4,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	80 62 57
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	300 520 420	the last opposed some some some some some som	18.0 18.0 18.0	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 4.00	46 66 39 P
1,2,3,7,8,9-HxCDF Total HxCDF	200 10000	يمين جيني شواد الماد ال	18.0 18.0	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	39 110 68 1400		18.0 18.0 18.0 18.0	2,3,7,8-TCDD-37Cl4	0.20	92
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	640 99 1600	*****	18.0 18.0 18.0	Total 2,3,7,8-TCDD Equivalence: 1800 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	490 1000		18.0 18.0			
OCDF OCDD	210 2900	****	36.0 36.0			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-07 Lab Sample ID 1034231004 Filename F60629A 11 Injected By BAL **Total Amount Extracted** 20.9 g Matrix Solid % Moisture 74.0 Dilution NA Dry Weight Extracted 5.42 g Collected 06/21/2006 **ICAL Date** 05/31/2006 Received 06/23/2006 CCal Filename(s) F60629A_01 & F60629A 16 Extracted 06/27/2006 Method Blank ID BLANK-10052 Analyzed 06/29/2006 22:32

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	32 3700	445 MM 445 445 445	11.0 A 0.37	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	82 91 68
2,3,7,8-TCDD Total TCDD	12 270	44 wheel each each	2.0 A 0.37	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00	80 101 123
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	160 14000	140	4.4 EA 4.2 A 1.8	1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00	83 84 90
1,2,3,7,8-PeCDD Total PeCDD	41 530	100 AN 100 AN 100	1.8 1.8	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	96 75 56 43
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	210 290 230	77 TO 10 TO 10	2.9 A 6.1 A 3.6 A	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 4.00	62 34 P
1,2,3,7,8,9-HxCDF Total HxCDF	140 8300	ada dali ilibi ada ada	3.1 A 1.8	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	18 70 31 900		2.1 A 2.4 A 1.8 1.8	2,3,7,8-TCDD-37Cl4	0.20	115
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	350 69 870	uph dish jare peperap dish dish mili dap ran mili sare dan dap ran	1.8 3.4 A 1.8	Total 2,3,7,8-TCDD Equivalence: 220 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	180 400	**************************************	1.8 1.8			
OCDF OCDD	170 840		3.7 3.7			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-09
Lab Sample ID 1034231005
Filename F60629A_12
Injected By BAL
Total Amount Extracted 18.3 g

Total Amount Extracted 18.3 g Matrix
% Moisture 25.0 Dilution
Dry Weight Extracted 13.7 g Collecte

 Dry Weight Extracted
 13.7 g
 Collected
 06/21/2006

 ICAL Date
 05/31/2006
 Received
 06/23/2006

 CCal Filename(s)
 F60629A_01 & F60629A_16
 Extracted
 06/27/2006

 Method Blank ID
 BLANK-10052
 Analyzed
 06/29/2006
 23:22

				,,			
Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery	
2,3,7,8-TCDF Total TCDF	6.00 160.00		0.61 A 0.15	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00	86 88	
2,3,7,8-TCDD Total TCDD	0.52 12.00	alle log yan gaw gan	0.38 JA 0.15	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	69 82 106	
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	28.00 500.00	120 	1.00 EA 0.73 0.73	1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00	84 57 68 69	
1,2,3,7,8-PeCDD Total PeCDD	2.20 25.00		0.73 J 0.73	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	76 61 56 42	
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	11.00 12.00 12.00		0.73 0.73 0.73	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 4.00	61 37 P	
1,2,3,7,8,9-HxCDF Total HxCDF	5.30 320.00	الله على الله الله الله الله الله الله الله ال	0.73 0.73	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA	
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	1.80 4.80 2.80 51.00		0.73 J 0.73 0.73 J 0.73	2,3,7,8-TCDD-37Cl4	0.20	104	
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	40.00 5.10 120.00	mer also has had had. Mer also had had. Also had had had.	0.73 0.73 0.73	Total 2,3,7,8-TCDD Equivalence: 23 ng/Kg (Using ITE Factors)			
1,2,3,4,6,7,8-HpCDD Total HpCDD	74.00 170.00		0.73 0.73				
OCDF OCDD	80.00 810.00	Mad Adds and Made Made Adjustes under supp	1.50 1.50				

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

! = Interference

Solid

NA

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



Percent

Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-12 Lab Sample ID 1034231006 Filename U60701A 21 Injected By BAL

Total Amount Extracted 18.0 g Matrix Solid % Moisture 53.9 Dilution NA

Dry Weight Extracted 8.30 g Collected 06/22/2006 07/01/2006 ICAL Date Received 06/23/2006 CCal Filename(s) U60701A_10 & U60701A 23 Extracted 06/28/2006 Method Blank ID **BLANK-10057** Analyzed 07/02/2006 06:23

EMPC Native Conc LRL Internal ng's

Isomers	ng/Kg	ng/Kg	ng/Kg	Standards	Added	Recovery
2,3,7,8-TCDF Total TCDF	9.3 1400.0	407 MG MA MA	8.0 A 0.24	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C	2.00 2.00	98 88
2,3,7,8-TCDD Total TCDD	7.3 160.0		4.5 A 0.24	1,2,3,7,8-PeCDF-13C 2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00	69 73 76
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	140.0 6900.0	170 	8.4 EA 7.3 A 1.2	1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00 2.00	142 P 116 105 105
1,2,3,7,8-PeCDD Total PeCDD	30.0 390.0		6.8 A 1.2	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	105 97 73 65
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	140.0 140.0 120.0 52.0		3.0 A 4.1 A 4.7 A	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 4.00	77 57
1,2,3,7,8,9-HxCDF Total HxCDF	4600.0		4.2 A 1.2	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	15.0 35.0 16.0 470.0		5.6 A 4.2 A 4.0 A 1.2	2,3,7,8-TCDD-37Cl4	0.20	96
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	200.0 41.0 500.0		3.4 A 5.0 A 1.2	Total 2,3,7,8-TCDD Equivalence: 150 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	90.0 190.0		2.9 A 1.2			
OCDF OCDD	82.0 240.0	100 AM 100 AM 100	6.8 A 10.0 A			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2.3.7.8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-14
Lab Sample ID 1034231007
Filename P60630A_09
Injected By BAL
Total Amount Extracted 20.3 g

Total Amount Extracted 20.3 g
% Moisture 27.7

Dry Weight Extracted 14.7 g

ICAL Date 05/20/2006
CCal Filename(s) P60630A_02 & P60630A_17
Method Blank ID BLANK-10057

 0.3 g
 Matrix
 Solid

 7.7
 Dilution
 NA

 4.7 g
 Collected
 06/22/2006

 5/20/2006
 Received
 06/23/2006

 60630A_02 & P60630A_17
 Extracted
 06/28/2006

Method Blank ID BLANK-10057 Analyzed 07/01/2006 01:23

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	600.0	4.9	0.450 EA 0.140	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00	89 89
2,3,7,8-TCDD Total TCDD	1.1 23.0	IPM THE MAN AND AND AND AND AND AND AND AND AND A	0.170 A 0.140	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00	67 76 98
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	8.1 190.0 2100.0	*** *** *** *** *** *** *** *** *** ***	0.920 A 0.870 A 0.680	1,2,3,4,7,6-HXCDF-13C 1,2,3,6,7,8-HXCDF-13C 2,3,4,6,7,8-HXCDF-13C 1,2,3,7,8,9-HXCDF-13C 1,2,3,4,7,8-HXCDD-13C	2.00 2.00 2.00 2.00 2.00	97 81 80 82
1,2,3,7,8-PeCDD Total PeCDD	3.3 41.0	MA AND MA AND MA	0.680 J 0.680	1,2,3,4,7,8-HXCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	94 76 79
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	17.0 30.0 75.0		0.680 0.680 0.680	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 2.00 4.00	75 106 80
1,2,3,7,8,9-HxCDF Total HxCDF	13.0 900.0		0.680 0.680	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	2.2 5.7 3.3 70.0	and this has not one. The black has been and one. The black has been and one. The black has been and one.	0.680 J 0.680 0.680 J 0.680	2,3,7,8-TCDD-37Cl4	0.20	95
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	35.0 6.5 87.0		0.680 0.680 0.680	Total 2,3,7,8-TCDD Equivalence: 110 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	21.0 44.0	We ald am like the	0.680 0.680			
OCDF OCDD	8.0 70.0	MA 140 MF 160 -160	1.400 1.400			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal


ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



Pace Analytical

Tel: 612-607-1700 Fax: 612- 607-6444

Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-16 Lab Sample ID 1034231008 Filename F60708A 05 Injected By BAL

Total Amount Extracted 16.1 g Matrix Solid % Moisture 17.0 Dilution NA

Dry Weight Extracted 13.3 g Collected 06/22/2006 ICAL Date 05/31/2006 Received 06/23/2006 F60708A_01 & F60708A 17 CCal Filename(s) Extracted 06/28/2006

Method Blank ID BLANK-10057 Analyzed 07/09/2006 00:09

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	ND ND	FF FF 600 AV. GA	0.150 0.150	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	92 88 74
2,3,7,8-TCDD Total TCDD	ND ND	ore and not the	0.150 0.150	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	85 103 89
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	ND ND ND		0.750 0.750 0.750	1,2,3,6,7,8-HXCDF-13C 2,3,4,6,7,8-HXCDF-13C 1,2,3,7,8,9-HXCDF-13C 1,2,3,4,7,8-HXCDD-13C	2.00 2.00 2.00 2.00 2.00	82 87 89 87
1,2,3,7,8-PeCDD Total PeCDD	ND ND	and and the same	0.750 0.750	1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	82 80 66
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF	ND ND ND ND		0.750 0.750 0.750 0.750	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C 1,2,3,4-TCDD-13C	2.00 2.00 4.00	82 64 NA
Total HxCDF	ND		0.750	1,2,3,7,8,9-HxCDD-13C	2.00	NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	ND ND ND ND	~~~~	0.750 0.750 0.750 0.750	2,3,7,8-TCDD-37Cl4	0.20	95
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	ND ND ND	70 VI 00 VI 100 WI 00 VI 00 VI WI 100 VI 00 VI	0.750 0.750 0.750	Total 2,3,7,8-TCDD Equivalence: 0.0044 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	ND ND	164 AND THE AND AND	0.750 0.750			
OCDF OCDD	ND 4.4	and the saw that das.	1.500 1.500 BJ			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal


ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



^Bace Analytical

Tel: 612-607-1700 Fax: 612-607-6444

Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-18 Lab Sample ID 1034231009 Filename P60630A 11 Injected By BAL

20.2 g **Total Amount Extracted** Matrix Solid % Moisture 82.7 Dilution NA

Dry Weight Extracted 3.50 g Collected 06/22/2006 ICAL Date 05/20/2006 Received 06/23/2006 CCal Filename(s) P60630A 02 & P60630A 17 Extracted 06/28/2006 Method Blank ID BLANK-10057 Analyzed 07/01/2006 02:59

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	1600.0	52	1.80 EA 0.57	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	82 77
2,3,7,8-TCDD Total TCDD	4.2 89.0	********	0.74 A 0.57	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	63 68 86
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	32.0 430.0 5500.0		4.50 A 3.20 A 2.90	1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00	92 70 72 74
1,2,3,7,8-PeCDD Total PeCDD	11.0 140.0	M M M M M	2.90 J 2.90	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	91 70 80 74
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	36.0 86.0 80.0	470 MP AND	5.00 A 3.00 A 2.90	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 4.00	109 76
1,2,3,7,8,9-HxCDF Total HxCDF	30.0 2300.0		3.20 A 2.90	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	9.5 25.0 17.0 300.0		2.90 J 2.90 2.90 2.90	2,3,7,8-TCDD-37Cl4	0.20	83
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	140.0 18.0 330.0	44 Me der ere per	2.90 2.90 2.90	Total 2,3,7,8-TCDD Equivalence: 260 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	280.0 630.0	100 and 100 and 100	2.90 2.90			
OCDF OCDD	87.0 1800.0	de fer en se se de se se se	5.70 5.70			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No....1034231



Method 8290 Analysis Results

Client - ESS Laboratory

Client's Sample ID 0606373-20 Lab Sample ID 1034231010 Filename P60630A 12 Injected By BAL **Total Amount Extracted** 16.0 g Matrix Solid % Moisture 12.3 Dilution NA Dry Weight Extracted 14.0 g Collected 06/22/2006 **ICAL** Date 05/20/2006 Received 06/23/2006 CCal Filename(s) P60630A 02 & P60630A 17 Extracted 06/28/2006 Method Blank ID BLANK-10057 Analyzed 07/01/2006 03:47

Native Isomers	Conc ng/Kg	EMPC ng/Kg	LRL ng/Kg	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF Total TCDF	0.24 3.70	dade dans mile bade bade bade dade dade dade dade dad	0.16 JA 0.14	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00	84 83
2,3,7,8-TCDD Total TCDD	ND 0.31	460 100 ma tạo dan	0.14 0.14 J	1,2,3,7,8-PeCDF-13C 2,3,4,7,8-PeCDD-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00	70 76 98
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF Total PeCDF	ND 0.73 5.60	ALL ME HOU MAY.	0.71 0.71 J 0.71	1,2,3,4,7,8-HXCDF-13C 1,2,3,6,7,8-HXCDF-13C 2,3,4,6,7,8-HXCDF-13C 1,2,3,7,8,9-HXCDF-13C 1,2,3,4,7,8-HXCDD-13C	2.00 2.00 2.00 2.00 2.00	97 78 79 82
1,2,3,7,8-PeCDD Total PeCDD	ND ND	** ** ** ***	0.71 0.71	1,2,3,4,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	98 73 78
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF	ND ND ND ND		0.71 0.71 0.71 0.71	1,2,3,4,7,6,9-npodr-13C 1,2,3,4,6,7,8-HpCDD-13C OCDD-13C 1,2,3,4-TCDD-13C	2.00 2.00 4.00	76 110 82 NA
Total HxCDF	3.90	~~~	0.71	1,2,3,7,8,9-HxCDD-13C	2.00	NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD Total HxCDD	ND ND ND 1.10		0.71 0.71 0.71 0.71 J	2,3,7,8-TCDD-37Cl4	0.20	85
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF Total HpCDF	2.00 ND 2.00	and the line and sub- and the line and spy	0.71 J 0.71 0.71 J	Total 2,3,7,8-TCDD Equivalence: 0.55 ng/Kg (Using ITE Factors)		
1,2,3,4,6,7,8-HpCDD Total HpCDD	7.40 13.00	की के का का का की का का का	0.71 0.71			
OCDF OCDD	3.10 64.00	*** *** *** **** **** **** ****	1.40 J 1.40			

Results reported on a dry weight basis

Conc = Concentration (Totals include 2,3,7,8-substituted isomers)

EMPC = Estimated Maximum Possible Concentration

A = Detection Limit based on signal-to-noise measurement

Nn = Value obtained from additional analysis

LRL = Lower Reporting Limit

I = Interference

E = PCDE Interference

S = Saturated signal

ND = Not Detected

NA = Not Applicable

NC = Not Calculated

* = See Discussion

Report No.....1034231

J = Concentration detected is below the calibration range

B = Less than 10 times higher than method blank level

P = Recovery outside of target range



Method 8290 Spike Sample Results

Client - ESS Laboratory

Client's Sample ID

Lab Sample ID Filename

Total Amount Extracted **ICAL** Date

CCal Filename(s) Method Blank ID

0606373-20 1034231010-MS

P60630A 14 16.2 g

05/20/2006

P60630A_02 & P60630A_17 BLANK-10057

Matrix Dilution Extracted Solid NA

06/28/2006 Analyzed 07/01/2006 05:23

Injected By BAL

Native Isomers	Qs (ng)	Qm (ng)	% Rec.	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF	0.20	0.18	92	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	73 75
2,3,7,8-TCDD	0.20	0.19	97	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	69 75 98
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF	1.00 1.00	1.11 1.06	111 106	1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00 2.00	81 65 66 70
1,2,3,7,8-PeCDD	1.00	0.92	92	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	78 65 65 59
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	1.00 1.00 1.00	1.00 1.09 1.08	100 109 108	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 2.00 4.00	87 53
1,2,3,7,8,9-HxCDF	1.00	1.07	107	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD	1.00 1.00 1.00	1.06 1.07 1.14	106 107 114	2,3,7,8-TCDD-37Cl4	0.20	79
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF	1.00 1.00	1.13 1.12	113 112			
1,2,3,4,6,7,8-HpCDD	1.00	0.98	98			
OCDF OCDD	2.00 2.00	1.70 2.18	85 109			

Qs = Quantity Spiked

Qm = Quantity Measured

Rec. = Recovery (Expressed as Percent)

P = Recovery outside of target range of 40-135%

X = Background subtracted value

E = PCDE Interference

Nn = Value obtained from additional analysis

NA = Not Applicable

* = See Discussion

Report No....1034231



Method 8290 Spike Sample Results

Client - ESS Laboratory

Client's Sample ID

Lab Sample ID Filename

Total Amount Extracted ICAL Date

CCal Filename(s)
Method Blank ID

0606373-20

1034231010-MSD P60630A 15

16.0 g 05/20/2006

05/20/2006 P60630A_02 & P60630A_17 BLANK-10057 Matrix Dilution Extracted Solid NA

Extracted 06/28/2006 Analyzed 07/01/2006 06:11

Injected By BAL

Native Isomers	Qs (ng)	Qm (ng)	% Rec.	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF	0.20	0.17	87	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00	75 78
2,3,7,8-TCDD	0.20	0.19	96	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C	2.00 2.00 2.00	64 76 101
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF	1.00 1.00	1.04 1.00	104 100	1,2,3,4,7,8-HxCDF-13C 1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C 1,2,3,4,7,8-HxCDD-13C	2.00 2.00 2.00 2.00	81 69 69 74
1,2,3,7,8-PeCDD	1.00	0.86	86	1,2,3,4,7,6-HXCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00 2.00	86 73 74 68
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF	1.00 1.00 1.00 1.00	0.95 1.05 1.02 1.03	95 105 102 103	1,2,3,4,6,7,8-HpCDD-13C OCDD-13C 1,2,3,4-TCDD-13C	2.00 2.00 4.00	92 68
1,2,0,1,0,0-110001	1.00	1.05	103	1,2,3,7,8,9-HxCDD-13C	2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD	1.00 1.00 1.00	1.11 0.98 1.02	111 98 102	2,3,7,8-TCDD-37Cl4	0.20	81
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF	1.00 1.00	1.08 1.13	108 113			
1,2,3,4,6,7,8-HpCDD	1.00	0.94	94			
OCDF OCDD	2.00 2.00	1.79 2.05	90 102			

Qs = Quantity Spiked

Qm = Quantity Measured

Rec. = Recovery (Expressed as Percent)

P = Recovery outside of target range of 40-135%

X = Background subtracted value

E = PCDE Interference

Nn = Value obtained from additional analysis

NA = Not Applicable

* = See Discussion

Report No....1034231

Pace Analytical Services, Inc. 1700 Elm Street - Suite 200 Minneapolis, MN 55414 Tel: 612-607-1700 Fax: 612- 607-6444

Pace Analytical"

Method 8290 Spike Sample Results

Client - ESS Laboratory

Client Sample ID	0606373-20				<u>م</u>	Dry Weights	eiahts	
Lab Sample ID MS ID MSD ID	1034231010 1034231010-MS 1034231010-MSD	Sam MS F MSD	Sample Filename MS Filename MSD Filename	P60630A_12 P60630A_14 P60630A_15	0A_12 0A_14 0A_15	Sample Ame MS Amount MSD Amount	Sample Amount 14.0 g MS Amount 14.2 g MSD Amount 14.0 g	
	Sample Conc. ng/Kg	MS/MSD Qs (ng)	MS Qm (ng)	MSD Qm (ng)	RPD	Backgroi MS % Rec.	Background Subtracted Rec. MSD % Rec.	RPD
2,3,7,8-TCDF		0.20	0.18	0.17	6.4	06	86	5.0
300		0.20	0.19	0.19	1.2	26	96	1.5
PeCDF		1.00	1.1	1.04	6.9	17	104	0.0
PeCDF		1.00	1.06	1.00	6.1	105	66	6.2
PeCDD		1.00	0.92	0.86	8.9	92	98	6.8
8-TXCDF		1.00	1.00	0.95	4.8	100	95	4.8
8-HXCDF		1.00	1.09	1.05	4.5	109	105	4.5
8-HXCDF		1.00	1.08	1.02	5.4	108	102	5.4
9-HXCDF		1.00	1.07	1.03	4.3	107	103	4.3
3-HXCDD		1.00	1.06	-	4.3	106	111	6.4
3-TXCDD		1.00	1.07	0.98	9.7	107	98	9.7
9-HXCDD		1.00	1.14	1.02	11.0	411	102	11.0
,8-HpCDF		1.00	1.13	1.08	4.5	110	105	4.5
3,9-HpCDF		1.00	1.12	1.13	1.0	112	113	1.0
,2,3,4,6,7,8-HpCDD		1.00	0.98	0.94	3.6	87	84	3.9
	3.148	2.00	1.70	1.79	5,2	83	87	5,4
	64.372	2.00	2.18	2.05	6.2	63	22	8.6
Definitions MS = Matrix Spike MSD = Matrix Spike Duplicate		CDD = Chlorinated dibenzo-p-dioxin	dibenzo-p-dioxin				The control of the co	THE OWN TO COME AND THE STREET AND T
Om = Quantity Measured Os = Quantity Sniked		T = Tetra Pe = Penta						
% Rec. = Percent Recovery RPD = Relative Percent Difference	Prence	Hx = Hexa Hn = Henta						
		0 = Octa						



Method 8290 Laboratory Control Spike Results

Client - ESS Laboratory

Lab Sample ID Filename

Total Amount Extracted ICAL Date

CCal Filename(s) Method Blank ID

LCS-10053 U60629A 05

10.1 g 02/13/2006

U60629A 04 & U60629A 20 BLANK-10052

Matrix Dilution Extracted Solid NA

06/27/2006 Analyzed 06/29/2006 12:56

Injected By SMT

Native Isomers	Qs (ng)	Qm (ng)	% Rec.	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF	0.20	0.21	105	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C	2.00 2.00	61 78
2,3,7,8-TCDD	0.20	0.17	85	1,2,3,7,8-PeCDF-13C 2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C	2.00 2.00 2.00	59 64 71
1,2,3,7,8-PeCDF 2,3,4,7,8-PeCDF	1.00 1.00	1.08 1.02	108 102	1,2,3,4,7,8-HxCDF-13C 1,2,3,6,7,8-HxCDF-13C 2,3,4,6,7,8-HxCDF-13C	2.00 2.00 2.00	85 70 74
1,2,3,7,8-PeCDD	1.00	0.88	88	1,2,3,7,8,9-HxCDF-13C 1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C	2.00 2.00 2.00 2.00	80 74 72 71
1,2,3,4,7,8-HxCDF 1,2,3,6,7,8-HxCDF 2,3,4,6,7,8-HxCDF	1.00 1.00 1.00	0.86 0.95 0.90	86 95 90	1,2,3,4,7,8,9-HpCDF-13C 1,2,3,4,6,7,8-HpCDD-13C OCDD-13C	2.00 2.00 4.00	64 81 64
1,2,3,7,8,9-HxCDF	1.00	0.91	91	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD	1.00 1.00 1.00	1.07 1.03 1.06	107 103 106	2,3,7,8-TCDD-37CI4	0.20	78
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF	1.00 1.00	1.00 0.99	100 99			
1,2,3,4,6,7,8-HpCDD	1.00	0.86	86			
OCDF OCDD	2.00 2.00	1.91 1.92	96 96			

Qs = Quantity Spiked Qm = Quantity Measured


Rec. = Recovery (Expressed as Percent) P = Recovery outside of target range

X = Background subtracted value

Nn = Value obtained from additional analysis

NA = Not Applicable * = See Discussion

Report No....1034231



Method 8290 Laboratory Control Spike Results

Client - ESS Laboratory

Lab Sample ID Filename **Total Amount Extracted**

ICAL Date CCal Filename(s) Method Blank ID

<u> Pace Analytical</u>

LCS-10058 F60630A 02 10.0 g 05/31/2006

F60630A_01 & F60630A_13 BLANK-10057

Matrix Dilution Extracted Analyzed

Solid NA

06/28/2006 06/30/2006 09:54

Injected	Ву	CVS
mjecteu	Бу	CVS

Native Isomers	Qs (ng)	Qm (ng)	% Rec.	Internal Standards	ng's Added	Percent Recovery
2,3,7,8-TCDF	0.20	0.17	86	2,3,7,8-TCDF-13C 2,3,7,8-TCDD-13C 1,2,3,7,8-PeCDF-13C	2.00 2.00 2.00	70 69
2,3,7,8-TCDD	0.20	0.18	90	2,3,4,7,8-PeCDF-13C 1,2,3,7,8-PeCDD-13C 1,2,3,4,7,8-HxCDF-13C	2.00 2.00 2.00 2.00	65 71 96
1,2,3,7,8-PeCDF	1.00	0.99	99	1,2,3,6,7,8-HxCDF-13C	2.00	78 73
2,3,4,7,8-PeCDF	1.00	0.91	91	2,3,4,6,7,8-HxCDF-13C 1,2,3,7,8,9-HxCDF-13C	2.00 2.00	79 76
1,2,3,7,8-PeCDD	1.00	0.81	81	1,2,3,4,7,8-HxCDD-13C 1,2,3,6,7,8-HxCDD-13C 1,2,3,4,6,7,8-HpCDF-13C 1,2,3,4,7,8,9-HpCDF-13C	2.00 2.00 2.00	83 76 72
1,2,3,4,7,8-HxCDF	1.00	0.87	87	1,2,3,4,6,7,8-HpCDD-13C	2.00 2.00	61 82
1,2,3,6,7,8-HxCDF	1.00	0.91	91	OCDD-13C	4.00	62
2,3,4,6,7,8-HxCDF 1,2,3,7,8,9-HxCDF	1.00 1.00	0.92 0.91	92 91	1,2,3,4-TCDD-13C 1,2,3,7,8,9-HxCDD-13C	2.00 2.00	NA NA
1,2,3,4,7,8-HxCDD 1,2,3,6,7,8-HxCDD 1,2,3,7,8,9-HxCDD	1.00 1.00 1.00	0.96 1.00 1.01	96 100 101	2,3,7,8-TCDD-37Cl4	0.20	68
1,2,3,4,6,7,8-HpCDF 1,2,3,4,7,8,9-HpCDF	1.00 1.00	1.00 1.05	100 105			
1,2,3,4,6,7,8-HpCDD	1.00	0.84	84			
OCDF OCDD	2.00 2.00	1.71 1.77	86 88			

Qs = Quantity Spiked Qm = Quantity Measured

Rec. = Recovery (Expressed as Percent) P = Recovery outside of target range X = Background subtracted value

Nn = Value obtained from additional analysis

NA = Not Applicable * = See Discussion

Report No....1034231

CHAIN OF CUSTODY

Page 3 of 3

ESS Laboratory

	€ (Y	1			-	5.1		4/5	70t.		<u> </u> ><	× 	×	×	×	×	<u> </u> ×	× ×		F-Filters
	25 LAW PROJECT ID 1606. 273	\ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \	Other	ĺ	lysis sisti	151243	w/a:	4/16 77	, v., .	×	×	×	X	×	×	×		X	V	W-Wipes F-F
0 0		Yes	PDF		Required Ana		MECT WBCT	W	NCLP-RCR	×	×	×	*	×	×	×	×	×	\forall	VIIO-0
	g raints	Electronic Deliverable	Format: Excel Access		Circle and/or Write Required Analysis	H	28 Vd 809	8082 6082 8082	FCRAS 8270 8270 8081 Posticides				****							DW-Drinking Water
10000	Tepotting Little	Electron	Format		5 - -	HeE	1 090 \$108 " 0: \$1	EbH 8100 EEX GB	Ebh											SW-Surface Water
÷	s required #	Other	Other				şış	nismoC	Code Number	1/62 (51/4	5.	# 6 "."	, r	12 12 12 12 12 12 12 12 12 12 12 12 12 1	5. "	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	e u	1, 2007		GW-Ground Water SW
Turn Time Standard Orher	er than 5 days, prior approval by labora where camples were collected from:	MA (RI CT NH NJ NY ME	Is this project for any of the following: MA-MCP Navy USACE	Project Name (20 Char. or le	Herham Sile		Zip PO#	Email Address	Sample Identification (20 Char or less)	SED 3101	Seb 3104	SED 30 01	SED 30 04	Sed 28 01	SED 28 03	SED 27 01	SED2904	SED 3201	SED 3204	Johd D-Studge WW-Waste Water
Inc.	185 Frances Avenue, Cranston, RI 02910-2211 Tel. (401) 461-7181 Fax (401) 461-4486	001-10		Project #	·	Address	**************************************	#	COMP COMP	X X	χS	X S	X S	X S	X 5	Х 5	χ 2	К 5	X 5	V-VOA Matrix: S-Soil SI
ngineering.	Cranston, RI 02910. Fax (401) 461.4486	, (101) m m	4. E.S. Programme of the Control of			el!	State	Fax #	Collection Time	1020	1050	1330	1345	1402	14.5	1442	1455	1522	1534	S-Sterile
of Thielsch E	185 Frances Avenue, (Tel. (401) 461-7181		or transport		下方がの	chris Ricardi		of 275 540	Date	6.21.06	ĮĮ.	li li	EF.	u u	a de la companya de l	H .	ж	*	11	pe: P-Poly G-Glass
Division	185 Franc Tel. (401)	www.essl		Co. Name	<u> </u>	Contact Person	Ŏ ŗ ġ	Telephone #	ESS LAB Sample#	_/	7,	γ-	7	5	3	N.	20	0	2	Container Type:

Please fax all changes to Chain of Custody in writing.

1 (White) Lab Copy 2 (Yellow) Client Receipt

Date/Time

Received by: (Signature)

Date/Time

Relinquished by: (Signature)

Date/Time

Received by: (Signature)

Date/Time

Relinquished by: (Signature)

4.8/10/8/14 Date/Time

(ed by (Signature)

4.21.06 RZ4 Date/Time

May Man Relinquished by: (Signature)

Relinquished by: (Signature)

Date/Time

Perceived by: (Signature)

4

Sampled by:

[] Pickup

No NA:

es

Seals Intact

Cooler Temp: 5.5

Comments:

[] Technicians

^{*}By circling MA-MCP, client acknowledges samples were collected in accordance with MADEP CAM VII A

		**************************************	<i>\$</i>	بر ۱۹۳۰ ا ک	(AT)	SAINOI Q	<u>ک</u>	-	×	molecularies	×	on the contraction	×	and transmitted the	×	***************************************	Djenstraanje E	anguarteasus	narayanasını	was property			turitia.
4 B M			····		7.7	10 C	· ×	<u> </u>		<u> </u>		×	-	_ >	1_	<u> </u>	ers	A COLUMN TO SERVICE			To the control of		
age 1 of 4 ESS LAB PROJECT ID CO 6373 No No #	4	5 repa	HW I	GI-71	1157				<u> </u>		 	- 	×		>	1	F-Filters	C. Miller of the Control of the Cont			Time	Lime	
STAB PROJ	11'	-	1.1	E-7 7 7	1	,	↓ ×		<u> </u>	 	<u> </u>	<u> </u>	×		×		ē				Date/Time	Date/Time	
	Ouner nalysis	.				1 AOV	<u> ×</u>		<u></u> ≥	<u> </u>	×	×		×			W-Wipes				12	14	
Page ESST Yes	d An	<u> </u>				AOV	×	×	×	×	×	X	×	×	×	×	ē.	9				+	
	Tire L	(E)	1) STV.	H/M LEVI-d:	ЭW	MELVEZ (13) MCb"											O-Oil	ZnAct,			rre)	ire)	
	Re l'	ļ		4DEN	8	TCLP-RCRA											i	8- Zr			gnati	gman	
ole A crosses	Circle and/or Write Required Analysis	F23		PP13	BCRA8	RCRAS	!										DW-Drinking Water	cid			Received by: (Signature)	Received by: (Signature)	
	d/oir		OLX H/V	⁄d 	979	voas 0/28											inkin	bic A			ved b	ved b	
	le an	8(8)	39 °	epicijse 809	9082 bcs b	Pesticides Pesticides											V-D	Asor			Recei	Recei	
Reporting Limits Electronic Delive	G. C.		Diesel BbH	*	\$HV4/M Hd3	HV4) O/M Hd9								·				£, 7-				-	
Repol Repol			<u> </u>	ONG 8013	100			 									Vater	leO!			ne	ne	
			id∧ Id∧	······		8021 WTBE/BT	1	1									ace \	6- N	4		Date/Time	Date/Time	
0		 				0888	 	 	-								SW-Surface Water	1OH,	1 4.5		Ω _a	Dat	
				1.2	enianno	Type of C	1/2//	-	<u></u>		<i>(,</i>	2	à	۵	w	4		5- N					
		***************************************		9131717	21102 10	Number	1	-							2	<u>_</u>	GW-Ground Water	Š.	1 .5		Relinquished by: (Signature)	Relinquished by: (Signature)	
r require	<u> </u>	—		T	T	Code	1+-	N.	1	3	14-	N	<u> </u>	~	Lf-	10	punc	4. H	1		(Sig)	(Sigr	
Other tory is	ĺ	-				Pres											%-Gr	SO,			d by:	d by:	
ME	3																- 1	3- H	9		uishe	uíshe	į
Light late		الب				(\$\$								ŀ			Water	£C,	2000		puilo	elinq	
andar roval cred I	4 0	N 지				ar. or ie								-			'aste'	2-]	2.		C	124	-
Turn Time Turn Time Standard If faster than 5 days, prior approval by laboratory is required # State where samples were collected from: MA RI OT NH NJ NY ME Other Electronic De Romar Product From Narry Change Collected from Narry Change Chan	Project Name (20 Clar. or less)	£ .			8	Sample Identification (20 Chac or less)							***************************************				WW-Waste Water	Preservation Code: 1- NP, 2- HCl, 3- H.SO., 4- HNO., 5- NaOH, 6- MeOH, 7- Asorbic Acid, 8-	Brigh		nne // //	me	NORTH BEAT
s were NH NA	Z -	Gorham		PO#	Email Address	tion						l					ŀ	ode:	150		Date/Time	Date/Time	_
mple St day	rojec	Ĕ			liai A	ntifica											D-Sludge	ion C	by:	įų.	Dave/	Δ̈́	
91909	<u> </u>				<u> </u>	I Ide		12		4	****	W	_	M		W	j	ervat	Sampled by:	iments:	3	····	
Turn Time If faster than State where ss MA KI Is this project MA-MCP					Î	amp	SEDZEOI	SEB 2003	SEDITO	SED IT 04	SED 1001	SEDIODS	6	Sed llez	SED 12 01	SED1203	SD-Solid	Pres	Sam	Comm			***************************************
Star Is I	4			Zip		S	52.1	2.0		Δ	Д	آم	SEDII	٦	٦	ā	S			1	y: (Signature)	ıture)	
	Project #	Address					NE L	N N	W.	\ <u>\</u>	8	5	2	N	N.	Š	S-Soil	4			y: (Signature	(Signe	l
7	Pro	- Ye									,						. S	Internal Use Only		[] Technicians_	279	deceived by: (Signature)	Manual Principles
0-2						MATRIX	悉		× 25.	8	**************************************	35	35.		3	35	Matrix:		ickuj	echni	Received	cervec	
291 448					1	CEVE	×	\times	\cong	X	<u>×</u>	×	\geq	×	×	ᅶ		Interi	[] Pickup	Ţ	Εğ.	<u> </u>	
ngineering, Inc. Cranston, RI 02910. Fax (401) 461-4486					#	COMB											V.VOA				4		I
1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1			اچــُ	State	Fax#	ion	0	6	,	0	M	1	1			- 1	ł		- 		Date/Time	£ime	
ee73 18to (40			A.	<i>V</i> 3		Collection Time	0840	5500	5160	0660	0945	0955	1015	1030	10 11	<u>=</u>	S-Sterile	No No	No NA:		ate/	Date/ Ime	-
Pax Pax	. 1		Kicardi		5401	0	೦	0	6	۵	0	J	****			1	•	z I	Z		Date/Time	-: -	
	MACTEC		·		1												G-Glass	λ. 	şs				
elsci emu 71817 ory.	7		Chris		芫	Date	0							İ			- N	Yes Fes	Yes		nunce)	(care)	
775 Av 61-,	Ž				7.97	ı	6.22	=	۵	2	=	٤	3	2		<u>:</u> ا	P-Poly	<u>ا</u> ف	1	\mathcal{U}	Sign.	ngic)	
of of the state of		- 1	-		1		ڻ									- 1	ن	ī	(4	- Y	· śc	
Gray (401 Rick)	ame	at Per	-		one #	AB		T	T	T	T						ner 1)	rese	tact	lemp		121161	THE COLUMN
Division of Thielsch Engineering, 185 Frances Avenue, Cranston, Tel. (401) 461-7181 Fax (401) www.esslaboratory.com	Co. Name	Contact Person		City	Telephone #	ESS LAB Sample#											Container Type:	Cooler (Tesent	Seals Intact	Cooler Temp: 2.	Relinquished by: (Signature)	<u> </u>	
		ΙÚ		<u> </u>	-		2	2	7	<u>~</u>	5 1		20			() [) J	S S	ပိ	e 1] e		a a